

Product datasheet for **HP223388**

ELP6 Human qPCR Primer Pair (BC014204)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54859
Forward Sequence:	GTCGTCTTCCAGGCTCAAAGG
Reverse Sequence:	CTCTCCACTGTCTACTGGCTTC
ACCN:	BC014204 , BC014204.2 , BC000623 , BC000623.2
UniProt ID:	Q0PNE2
Synonyms:	chromosome 3 open reading frame 75; FLJ20211; TMEM103; transmembrane protein 103
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)