

Product datasheet for **HP223394**

CHMP4C Human qPCR Primer Pair (BC014321)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	92421
Forward Sequence:	AGACTGAGGAGATGCTGGGCAA
Reverse Sequence:	TAGTGCCTGTAATGCAGCTCGC
ACCN:	BC014321 , BC014321.1 , BQ778244 , BX492886
UniProt ID:	Q96CF2
Synonyms:	chromatin modifying protein 4C; MGC108776; Shax3; Snf7 homologue associated with Alix 3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)