

Product datasheet for **HP223433**

PINX1 Human qPCR Primer Pair (BC015479)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54984
Forward Sequence:	TCCAGAGGAGAACGAAACCACG
Reverse Sequence:	TAGTTCCACCTGCGTCTCAGA
ACCN:	<u>BC015479</u> , <u>BC015479.1</u> , <u>BC093762</u> , <u>BC143382</u>
UniProt ID:	<u>Q96BK5</u>
Synonyms:	67-11-3 protein; FLJ20565; hepatocellular carcinoma-related putative tumor suppressor; LPTL; LPTS; MGC8850; PIN2-interacting protein 1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)