

Product datasheet for **HP223512**

GSTM2 Human qPCR Primer Pair (BC017836)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2946
Forward Sequence:	AGATCACCCAGAGCAACGCCAT
Reverse Sequence:	GGCTGTCCATAAACTGGTTCTCC
ACCN:	<u>BC017836</u> , <u>BC110380</u> , <u>BC110380.1</u> , <u>BC105038</u> , <u>BC105066</u> , <u>BM924181</u> , <u>BT019947</u>
UniProt ID:	<u>P28161</u>
Synonyms:	GST4; GSTM; GSTM2-2; GTHMUS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)