

Product datasheet for **HP223553**

C8orf55 (THEM6) Human qPCR Primer Pair (BC019104)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51337
Forward Sequence:	GTACCTGGTGCGCCTTCCGTG
Reverse Sequence:	GCGCGTTGTTTCATGTGCAGCAG
ACCN:	BC019104 , BC001311 , BC001311.1 , BM744164
UniProt ID:	Q8WUY1
Synonyms:	chromosome 8 open reading frame 55; DSCD75; mesenchymal stem cell protein DSCD75
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)