

Product datasheet for **HP223573**

SNX12 Human qPCR Primer Pair (BC020559)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	29934
Forward Sequence:	CCGCCAAGTAACTTCCTGGAGA
Reverse Sequence:	GTACGCAGGACTCCTTTAGCTTG
ACCN:	BC020559 , BC033809 , BC103847 , BC024301 , BC044059 , BC063716 , BC080535 , BC103848 , BC108323 , BM805683 , BT007203
UniProt ID:	Q9UMY4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)