

Product datasheet for **HP223698**

WDR51B (POC1B) Human qPCR Primer Pair (BC026080)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	282809
Forward Sequence:	GAACCGTGAGACTCTGGATTCC
Reverse Sequence:	CGCTGGCGATACATGCTCCATA
ACCN:	BC026080 , BC026080.2 , BU521431 , BX648293
UniProt ID:	Q8TC44
Synonyms:	FLJ14923; FLJ41111; TUWD12; WD repeat domain 51B
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)