

## Product datasheet for **HP225306**

### JRK Human qPCR Primer Pair (NM\_003724)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8629
Forward Sequence:	GCATCCAGTGAAAAGCAGTCAGC
Reverse Sequence:	CATCAGCGTTGTAACCTGCTCG
ACCN:	<a href="#">NM_003724</a> , <a href="#">NM_003724.1</a> , <a href="#">NM_003724.2</a> , <a href="#">NM_003724.3</a> , <a href="#">BC043351</a> , <a href="#">BC024297</a> , <a href="#">BC030566</a> , <a href="#">BC036794</a> , <a href="#">BC071862</a> , <a href="#">BX640755</a> , <a href="#">NM_003724.4</a>
UniProt ID:	<a href="#">O75564</a>
Synonyms:	jerky; JH8
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)