

Product datasheet for **HP225314**

KLHL9 Human qPCR Primer Pair (NM_018847)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55958
Forward Sequence:	CCACTGATGACACCACAGGATC
Reverse Sequence:	CTATCTGACTGCATCACTGGCTG
ACCN:	NM_018847 , NM_018847.1 , NM_018847.2 , NM_018847.3 , BC039133 , BC039133.2 , BC052776 , BC064958 , BC105008 , BC113513 , BX538121 , NM_018847.4
UniProt ID:	Q9P2J3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)