

Product datasheet for **HP225348**

NAPE PLD (NAPEPLD) Human qPCR Primer Pair (NM_198990)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	222236
Forward Sequence:	GCAGGAGATACTGGTTATTGCC
Reverse Sequence:	GGACATCAGTGTGAATCCTTACAG
ACCN:	<u>NM_198990</u> , <u>NM_198990.1</u> , <u>NM_198990.2</u> , <u>NM_198990.3</u> , <u>NM_198990.4</u> , <u>BC071604</u> , <u>BC071604.1</u> , <u>BC037350</u> , <u>BC057294</u>
UniProt ID:	<u>Q6IQ20</u>
Synonyms:	C7orf18; FMP30; NAPE-PLD
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)