

## Product datasheet for **HP225390**

### PKIG Human qPCR Primer Pair (NM\_181805)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11142
Forward Sequence:	TGACATCCAGGGAGACTCAGAG
Reverse Sequence:	CTGGTTGCCAGCTTCCTTGCT
ACCN:	<u><a href="#">NM_181805</a></u> , <u><a href="#">NM_181805.1</a></u> , <u><a href="#">NM_181805.2</a></u> , <u><a href="#">BC104256</a></u> , <u><a href="#">BC104257</a></u> , <u><a href="#">BF569812</a></u> , <u><a href="#">BG388302</a></u> , <u><a href="#">BG717090</a></u> , <u><a href="#">BI461377</a></u> , <u><a href="#">BQ010165</a></u> , <u><a href="#">BU161358</a></u>
UniProt ID:	<u><a href="#">Q9Y2B9</a></u>
Synonyms:	PKI-gamma
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)