

Product datasheet for **HP225421**

ERGIC3 Human qPCR Primer Pair (NM_198398)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51614
Forward Sequence:	GCTATGGTGCTGAGGCAGAAGA
Reverse Sequence:	CTCCTGCATCTTCTGGCTGAAG
ACCN:	<u>NM_198398</u> , <u>NM_198398.1</u> , <u>BC014014</u> , <u>BC014014.2</u> , <u>BC009765</u> , <u>BU528333</u> , <u>NM_198398.2</u>
UniProt ID:	<u>Q9Y282</u>
Synonyms:	C2orf47; C20orf47; CGI-54; dj477O4.2; Erv46; NY-BR-84; PRO0989; SDBCAG84
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)