

Product datasheet for **HP225444**

CSNK1G3 Human qPCR Primer Pair (NM_004384)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1456
Forward Sequence:	CACAAGAGCCTTACAGGAACAGC
Reverse Sequence:	AAGCCTTGCCAAGGAAGACTGC
ACCN:	<u>NM_004384</u> , <u>NM_004384.1</u> , <u>NM_004384.2</u> , <u>NM_004384.3</u> , <u>NM_004384.4</u> , <u>BC047567</u> , <u>BC047567.1</u> , <u>BX116822</u>
UniProt ID:	<u>Q9Y6M4</u>
Synonyms:	CKI-gamma 3; CSNK1G3L
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)