

Product datasheet for **HP225450**

SLC35A2 Human qPCR Primer Pair (NM_005660)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7355
Forward Sequence:	ATGGCGGAAGTGCTCAAAGGTC
Reverse Sequence:	CTTGAGCGTGTCCACATACTGC
ACCN:	BC035747 , NM_005660 , NM_005660.1 , NM_005660.2 , BC035747.1 , BI820134 , NM_005660.3
UniProt ID:	P78381
Synonyms:	CDG2M; CDGX; UDP-Gal-Tr; UGALT; UGAT; UGT; UGT1; UGT2; UGTL
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)