

Product datasheet for **HP225459**

RNF220 Human qPCR Primer Pair (NM_018150)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55182
Forward Sequence:	ACCGATGACCTCCACCATTTCAG
Reverse Sequence:	CTGCCCTTCATCTTGCTTCTCTC
ACCN:	NM_018150 , NM_018150.1 , NM_018150.2 , NM_018150.3 , BC034221 , BC034221.1 , BC000279 , BC062600 , BC073835 , BC098266 , BC098300 , BC105790 , NM_018150.4
UniProt ID:	Q5VTB9
Synonyms:	C1orf164
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)