

Product datasheet for **HP225553**

MBNL2 Human qPCR Primer Pair (NM_144778)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10150
Forward Sequence:	ATTCTGCCCACCACGCCTGTTA
Reverse Sequence:	GGCACAGTTTCCTCGCTGGAAC
ACCN:	<u>NM_144778</u> , <u>NM_144778.1</u> , <u>NM_144778.2</u> , <u>NM_144778.3</u> , <u>BC104038</u> , <u>BC020418</u> , <u>BC020430</u> , <u>BC032426</u> , <u>BC104039</u> , <u>BC104040</u> , <u>BE246602</u> , <u>NM_144778.4</u>
UniProt ID:	<u>Q5VZF2</u>
Synonyms:	MBLL; MBLL39; PRO2032
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)