

## Product datasheet for **HP225569**

### POMZP3 Human qPCR Primer Pair (NM\_152992)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22932
Forward Sequence:	GCAAAGGAGACTGTACTGAGTGC
Reverse Sequence:	CAGAAGAGGCATCAGTGAGACC
ACCN:	<a href="#">BC017101</a> , <a href="#">NM_152992</a> , <a href="#">NM_152992.1</a> , <a href="#">NM_152992.2</a> , <a href="#">BC000487</a> , <a href="#">BC144677</a> , <a href="#">BC148538</a> , <a href="#">BC153130</a> , <a href="#">BG461659</a> , <a href="#">BI084456</a>
UniProt ID:	<a href="#">Q6PJE2</a>
Synonyms:	POM-ZP3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)