

Product datasheet for **HP225593**

ASRGL1 Human qPCR Primer Pair (NM_025080)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80150
Forward Sequence:	GCTGGAGGTTATGCCGACAATG
Reverse Sequence:	ATAGGTCCGCAGCCTCTTCTAC
ACCN:	<u>NM_025080</u> , <u>NM_025080.1</u> , <u>NM_025080.2</u> , <u>NM_025080.3</u> , <u>BC021295</u> , <u>BC021295.2</u> , <u>BC093070</u> , <u>BC006267</u> , <u>BC064963</u> , <u>BP308677</u> , <u>BX640832</u> , <u>NM_025080.4</u>
UniProt ID:	<u>Q7L266</u>
Synonyms:	ALP; ALP1; CRASH
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)