

## Product datasheet for **HP225606**

### ZXDC Human qPCR Primer Pair (NM\_025112)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79364
Forward Sequence:	AGGATGTGGGTGCTCCGAAAAG
Reverse Sequence:	CGGAGCTTCTAGCTGAGGTAAG
ACCN:	<u><a href="#">NM_025112</a></u> , <u><a href="#">NM_025112.1</a></u> , <u><a href="#">NM_025112.2</a></u> , <u><a href="#">NM_025112.3</a></u> , <u><a href="#">NM_025112.4</a></u> , <u><a href="#">BC002940</a></u> , <u><a href="#">BC012729</a></u> , <u><a href="#">BC063713</a></u> , <u><a href="#">BC073859</a></u> , <u><a href="#">BC110644</a></u> , <u><a href="#">BM686637</a></u> , <u><a href="#">NM_025112.5</a></u>
UniProt ID:	<u><a href="#">Q2QGD7</a></u>
Synonyms:	ZXDL
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)