

Product datasheet for **HP225614**

TP53BP2 Human qPCR Primer Pair (NM_005426)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7159
Forward Sequence:	AGCACTGGGAATGCTCTGGATC
Reverse Sequence:	GGCATTGGACTGGTCTACTGCA
ACCN:	<u>NM_005426</u> , <u>NM_005426.1</u> , <u>NM_005426.2</u> , <u>BC058918</u> , <u>BC058918.1</u> , <u>BC040247</u> , <u>BC046150</u> , <u>BX425756</u>
UniProt ID:	<u>Q13625</u>
Synonyms:	53BP2; ASPP2; BBP; P53BP2; PPP1R13A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)