

Product datasheet for **HP225621**

UHRF1 Human qPCR Primer Pair (NM_013282)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	29128
Forward Sequence:	GACAAGCAGCTCATGTGCGATG
Reverse Sequence:	AGTACCACCTCGCTGGCATCAT
ACCN:	NM_013282 , NM_013282.1 , NM_013282.2 , NM_013282.3 , NM_013282.4 , BC113875 , BC039136 , BC046179
UniProt ID:	Q96T88
Synonyms:	hNP95; hUHRF1; huNp95; ICBP90; Np95; RNF106; TDRD22
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)