

Product datasheet for **HP225623**

ENAH Human qPCR Primer Pair (NM_018212)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55740
Forward Sequence:	CGCCCTTAACTGGACTTGACAG
Reverse Sequence:	TAGATGAGCGGAGTTCACACC
ACCN:	NM_018212 , NM_018212.1 , NM_018212.2 , NM_018212.3 , NM_018212.4 , NM_018212.5 , BC095481 , BC095481.1 , BC010414 , BC037319 , BC065238 , BQ672102 , BU675547 , NM_018212.6
UniProt ID:	Q8N8S7
Synonyms:	ENA; MENA; NDPP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)