

Product datasheet for **HP225636**

GPRI77 (WLS) Human qPCR Primer Pair (NM_024911)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79971
Forward Sequence:	GTTGGCTCCTTCTGCCTCTTCA
Reverse Sequence:	AGGCAGATTCCAGCCACGATGA
ACCN:	<u>NM_024911</u> , <u>NM_024911.1</u> , <u>NM_024911.2</u> , <u>NM_024911.3</u> , <u>NM_024911.4</u> , <u>NM_024911.5</u> , <u>NM_024911.6</u> , <u>BC110826</u> , <u>BC007211</u> , <u>BC137109</u> , <u>BC137113</u> , <u>BG701224</u> , <u>BX537492</u> , <u>BX538320</u> , <u>BX648748</u> , <u>NM_024911.7</u>
UniProt ID:	<u>Q5T9L3</u>
Synonyms:	C1orf139; EVI; GPRI77; mig-14; MRP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

