

Product datasheet for **HP225643**

SLC30A5 Human qPCR Primer Pair (NM_024055)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64924
Forward Sequence:	TTTGAAGGCTGTGGGACTTTTCG
Reverse Sequence:	GGTGTTTGTAATAGTTTCCCAG
ACCN:	<u>NM_024055</u> , <u>NM_024055.2</u> , <u>NM_024055.3</u> , <u>NM_024055.4</u> , <u>BC000808</u> , <u>BC003411</u> , <u>BC008198</u> , <u>BC017441</u> , <u>BC130452</u> , <u>BC130454</u> , <u>BC144349</u> , <u>BC148650</u> , <u>BC153202</u> , <u>BX537394</u> , <u>NM_024055.5</u>
UniProt ID:	<u>Q8TAD4</u>
Synonyms:	ZnT-5; ZNT5; ZNTL1; ZTL1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)