

## Product datasheet for **HP225645**

### CRTC3 Human qPCR Primer Pair (NM\_022769)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64784
Forward Sequence:	CTCGATCCTGTGATGTTGGAGG
Reverse Sequence:	GTGTCGAGTAGTGGAGGTTGGT
ACCN:	<a href="#">NM_022769</a> , <a href="#">NM_022769.1</a> , <a href="#">NM_022769.2</a> , <a href="#">NM_022769.3</a> , <a href="#">NM_022769.4</a> , <a href="#">BC019060</a> , <a href="#">BC074730</a> , <a href="#">BC074731</a> , <a href="#">BP238198</a> , <a href="#">BP365179</a> , <a href="#">BQ945777</a> , <a href="#">NM_022769.5</a>
UniProt ID:	<a href="#">Q6UUU7</a>
Synonyms:	TORC-3; TORC3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)