

Product datasheet for **HP225653**

ZNF295 (ZBTB21) Human qPCR Primer Pair (NM_020727)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	49854
Forward Sequence:	CGGTCACAAGAGACCTGTCTCT
Reverse Sequence:	ACACAGGAGACCCGTGCTCATT
ACCN:	NM_020727 , NM_020727.1 , NM_020727.2 , NM_020727.3 , NM_020727.4 , BC063290 , BC063290.1 , BI767542 , BM893936 , BP230048
UniProt ID:	Q9ULJ3
Synonyms:	ZNF295
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)