

Product datasheet for **HP225698**

MOG1 (RANGRF) Human qPCR Primer Pair (NM_016492)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	29098
Forward Sequence:	GAGCCTGATAGTGGAACCTCTCG
Reverse Sequence:	TGAGAGGCTGAACAGACTCCAC
ACCN:	BC012552 , NM_016492 , NM_016492.1 , NM_016492.2 , NM_016492.3 , NM_016492.4 , BC012552.1 , BC100017 , BC006486 , NM_016492.5
UniProt ID:	Q9HD47
Synonyms:	HSPC165; HSPC236; MOG1; RANGNRF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)