

Product datasheet for **HP225722**

NIPP1 (PPP1R8) Human qPCR Primer Pair (NM_014110)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5511
Forward Sequence:	GAAGCCTCAGACATTGCCATCG
Reverse Sequence:	GGCAGTGTGAACTCTGTCAGG
ACCN:	<u>NM_014110</u> , <u>NM_014110.2</u> , <u>NM_014110.3</u> , <u>NM_014110.4</u> , <u>BC093017</u> , <u>BC093017.1</u> , <u>BC001597</u> , <u>BC003405</u> , <u>BC013360</u> , <u>BI459680</u> , <u>NM_014110.5</u>
UniProt ID:	<u>Q12972</u>
Synonyms:	ARD-1; ARD1; NIPP-1; NIPP1; PRO2047
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)