

Product datasheet for **HP225735**

SPFH2 (ERLIN2) Human qPCR Primer Pair (NM_007175)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11160
Forward Sequence:	TACGGGCAGAAGGTGATGGAGA
Reverse Sequence:	GCA GTGTAGCACTCAGCATCTG
ACCN:	<u>NM_007175</u> , <u>NM_007175.2</u> , <u>NM_007175.3</u> , <u>NM_007175.4</u> , <u>NM_007175.5</u> , <u>NM_007175.6</u> , <u>BC005950</u> , <u>BC035980</u> , <u>BC048308</u> , <u>BC050611</u> , <u>BC067765</u> , <u>BI560439</u> , <u>BP353279</u> , <u>BP363612</u> , <u>BU608317</u> , <u>NM_007175.8</u>
UniProt ID:	<u>O94905</u>
Synonyms:	C8orf2; Erlin-2; NET32; SPFH2; SPG18
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)