

Product datasheet for **HP225737**

Clathrin light chain (CLTA) Human qPCR Primer Pair (NM_007096)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1211
Forward Sequence:	CGATTGCAGTCAGAGCCTGAAAG
Reverse Sequence:	TAGCTGCTCGTCCTGTCTTGCA
ACCN:	NM_007096 , NM_007096.1 , NM_007096.2 , NM_007096.3 , BC009201 , BC009201.2 , BC019287 , BF797444 , BI090227 , BI600540 , BM423666 , BT007170 , NM_007096.4
UniProt ID:	P09496
Synonyms:	LCA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)