

## Product datasheet for **HP225766**

### **c Maf (MAF) Human qPCR Primer Pair (NM\_005360)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	4094
<b>Forward Sequence:</b>	AGAAGTTGGTGAGCAGCGGCTT
<b>Reverse Sequence:</b>	CACTGATGGCTCCAACTTGCGA
<b>ACCN:</b>	<u><a href="#">NM_005360</a></u> , <u><a href="#">NM_005360.2</a></u> , <u><a href="#">NM_005360.3</a></u> , <u><a href="#">NM_005360.4</a></u> , <u><a href="#">BC081542</a></u> , <u><a href="#">BC038438</a></u> , <u><a href="#">BE644986</a></u> , <u><a href="#">BF726289</a></u> , <u><a href="#">NM_005360.5</a></u>
<b>UniProt ID:</b>	<u><a href="#">O75444</a></u>
<b>Synonyms:</b>	AYGRP; c-MAF; CCA4; CTRCT21
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)