

Product datasheet for **HP225767**

PODXL Human qPCR Primer Pair (NM_005397)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5420
Forward Sequence:	CCTGAACCTCACAGGAAACACC
Reverse Sequence:	TGGAACAGATGCCAGCCGTATG
ACCN:	NM_005397 , NM_005397.1 , NM_005397.2 , NM_005397.3 , BC093730 , BC093730.1 , BC042466 , BC112035 , BC143318 , BG006197 , BP355358 , BX641124 , NM_005397.4
UniProt ID:	O00592
Synonyms:	gp135; Gp200; PC; PCLP; PCLP-1; PDX; PODXL1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)