

Product datasheet for **HP225781**

PLAA Human qPCR Primer Pair (NM_004253)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9373
Forward Sequence:	AGGCTGGAAGATGTGAGAGGAC
Reverse Sequence:	TCGCCAGTGATTTGCCACCTTC
ACCN:	<u>NM_004253</u> , <u>NM_004253.2</u> , <u>BC032551</u> , <u>BC032551.1</u>
UniProt ID:	<u>Q9Y263</u>
Synonyms:	DOA1; DOA1 homolog; DOA1, PLAP, PLA2P, FLJ11281, FLJ12699; FLJ11281; FLJ12699; OTTHUMP00000021164; OTTHUMP00000045176; phospholipase A2 activating protein; phospholipase A2-activating protein; PLA2P; PLAP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)