

Product datasheet for **HP225789**

INPP4B Human qPCR Primer Pair (NM_003866)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 8821 |
| Forward Sequence: | GATGCCATTCCCTCTCACCCAA |
| Reverse Sequence: | TTGGCTCGGTTCACTTCTC |
| ACCN: | NM_003866 , NM_003866.1 , NM_003866.2 , NM_003866.3 , BC110918 , BC110918.1 , BC005273 , BC005274 , BC017924 , BC072447 , BC133005 , BC144255 , BG708971 , BP360389 , BX649090 |
| UniProt ID: | O15327 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)