

Product datasheet for **HP225799**

ZNF207 Human qPCR Primer Pair (NM_003457)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7756
Forward Sequence:	ATGATGCCACCTGGACCAGGAA
Reverse Sequence:	GCTGAAACAGCCTGTGCTTGAG
ACCN:	<u>BC002372</u> , <u>BC008023</u> , <u>NM_003457</u> , <u>NM_003457.1</u> , <u>NM_003457.2</u> , <u>NM_003457.3</u> , <u>BC002372.2</u> , <u>BC008023.1</u> , <u>BC000962</u> , <u>BT007365</u>
UniProt ID:	<u>O43670</u>
Synonyms:	BuGZ; hBuGZ
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)