

## Product datasheet for **HP225806**

### Aurora C (AURKC) Human qPCR Primer Pair (NM\_003160)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6795
Forward Sequence:	CCACAGTGAGACTTACAGACGC
Reverse Sequence:	GCTGGTATCTGAGAAGCCTGGA
ACCN:	<u>NM_003160</u> , <u>NM_003160.1</u> , <u>NM_003160.2</u> , <u>BC075064</u> , <u>BX092186</u> , <u>NM_003160.3</u>
UniProt ID:	<u>Q9UQB9</u>
Synonyms:	AIE2; AIK3; ARK3; AurC; aurora-C; HEL-S-90; SPGF5; STK13
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)