

Product datasheet for **HP225811**

RAD1 Human qPCR Primer Pair (NM_002853)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5810
Forward Sequence:	GGATATGACGAGTGAAGTCCTAC
Reverse Sequence:	AGTCAAGGTGGAACTTCCTGC
ACCN:	<u>BC009804</u> , <u>NM_002853</u> , <u>NM_002853.1</u> , <u>NM_002853.2</u> , <u>NM_002853.3</u> , <u>BC009804.1</u> , <u>BC006837</u> , <u>BC006837.1</u> , <u>BC035771</u> , <u>BC037857</u> , <u>BT006908</u> , <u>NM_002853.4</u>
UniProt ID:	<u>O60671</u>
Synonyms:	HRAD1; REC1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)