

## Product datasheet for **HP225829**

### Estrogen Receptor beta (ESR2) Human qPCR Primer Pair (NM\_001437)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2100
Forward Sequence:	ATGGAGTCTGGTCGTGTGAAGG
Reverse Sequence:	TAACTTCCGAAGTCGGCAGG
ACCN:	<u>NM_001437</u> , <u>NM_001437.1</u> , <u>NM_001437.2</u> , <u>BC024181</u> , <u>BX094834</u> , <u>BX247957</u>
UniProt ID:	<u>Q92731</u>
Synonyms:	ER-BETA; Erb; ESR-BETA; ESRB; ESTRB; NR3A2; ODG8
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)