

## Product datasheet for **HP225832**

### CTBP1 Human qPCR Primer Pair (NM\_001328)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1487
Forward Sequence:	AGATGCCCATCCTGAAGGACGT
Reverse Sequence:	GAGGGCTTTGAACTTCTCCAGG
ACCN:	<a href="#">NM_001328</a> , <a href="#">NM_001328.1</a> , <a href="#">NM_001328.2</a> , <a href="#">BC053320</a> , <a href="#">BC098400</a> , <a href="#">BC011655</a> , <a href="#">BC064333</a> , <a href="#">BC072021</a> , <a href="#">BM989059</a> , <a href="#">BU553431</a> , <a href="#">NM_001328.3</a>
UniProt ID:	<a href="#">Q13363</a>
Synonyms:	BARS; HADDTS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)