

Product datasheet for **HP225834**

LARP1 Human qPCR Primer Pair (NM_033551)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23367
Forward Sequence:	TAGAGCGAGACTTCTTCCTGCG
Reverse Sequence:	ATCTCCACCACCTTGCTGTCCT
ACCN:	NM_033551 , NM_033551.2 , BC001460 , BC001460.2 , BC003063 , BC010144 , BC015573 , BC033856 , BC038220 , BQ003824 , BQ422667
UniProt ID:	Q6PKG0
Synonyms:	Lar1; LARP; Lhp1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)