

Product datasheet for **HP225852**

JNK1 (MAPK8) Human qPCR Primer Pair (NM_139049)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5599
Forward Sequence:	GACGCCTTATGTAGTGACTCGC
Reverse Sequence:	TCCTGGAAGAGGATTTTGTGGC
ACCN:	<u>NM_139049</u> , <u>NM_139049.1</u> , <u>NM_139049.2</u> , <u>NM_139049.3</u> , <u>BC130570</u> , <u>BC130572</u> , <u>BC144063</u> , <u>NM_139049.4</u>
UniProt ID:	<u>P45983</u>
Synonyms:	JNK; JNK-46; JNK1; JNK1A2; JNK21B1/2; PRKM8; SAPK1; SAPK1c
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)