

## Product datasheet for **HP225853**

### PPM1M Human qPCR Primer Pair (NM\_144641)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	132160
Forward Sequence:	CAGCTCAGAGTCCTGGACACAA
Reverse Sequence:	CCAGAGTCCATCAGTTGCCATG
ACCN:	<u><a href="#">NM_144641</a></u> , <u><a href="#">NM_144641.1</a></u> , <u><a href="#">NM_144641.2</a></u> , <u><a href="#">NM_144641.3</a></u> , <u><a href="#">BC009644</a></u> , <u><a href="#">BF516133</a></u>
UniProt ID:	<u><a href="#">Q96MI6</a></u>
Synonyms:	PP2C-eta; PP2CE; PP2Ceta
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)