

Product datasheet for **HP225854**

ADHFE1 Human qPCR Primer Pair (NM_144650)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	137872
Forward Sequence:	GTGAGAGTGGAACCAACGGATTC
Reverse Sequence:	AGCAGCCTTACAGGTGCCATG
ACCN:	NM_144650 , NM_144650.1 , NM_144650.2 , BC047492 , BC064634 , BC110290 , BE241883 , BG698892 , BQ182872 , NM_144650.3
UniProt ID:	Q8IWW8
Synonyms:	ADH8; HMFT2263; HOT
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)