

## Product datasheet for **HP225865**

### CD56 (NCAM1) Human qPCR Primer Pair (NM\_181351)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4684
Forward Sequence:	CATCACCTGGAGGACTTCTACC
Reverse Sequence:	CAGTGACTGGATGCTCTTCAGG
ACCN:	<a href="#">NM_181351</a> , <a href="#">NM_181351.1</a> , <a href="#">NM_181351.2</a> , <a href="#">NM_181351.3</a> , <a href="#">NM_181351.4</a> , <a href="#">BC047244</a> , <a href="#">BC047244.1</a> , <a href="#">BC014205</a> , <a href="#">BC019845</a> , <a href="#">BC029119</a>
UniProt ID:	<a href="#">P13591</a>
Synonyms:	CD56; MSK39; NCAM
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)