

Product datasheet for **HP225893**

TTYH1 Human qPCR Primer Pair (NM_020659)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57348
Forward Sequence:	GCTGTCTACCTCATCCGCTTCT
Reverse Sequence:	CTGTTGCCATAGAAACCGATGCC
ACCN:	BC019358 , NM_020659 , NM_020659.1 , NM_020659.2 , NM_020659.3 , BC011347 , BC011347.1 , BC065919 , BU737336 , BX402889 , NM_020659.4
UniProt ID:	Q9H313
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)