

## Product datasheet for **HP225917**

### **RUNX1 Human qPCR Primer Pair (NM\_001754)**

#### **Product data:**

|                          |   |
|--------------------------|---|
| <b>Product Type:</b>     | qPCR Primer Pairs   |
| <b>Gene ID:</b>          | 861   |
| <b>Forward Sequence:</b> | CCACCTACCACAGAGCCATCAA  |
| <b>Reverse Sequence:</b> | TTCACTGAGCCGCTCGGAAAAG  |
| <b>ACCN:</b>             | <u><a href="#">NM_001754</a></u> , <u><a href="#">NM_001754.1</a></u> , <u><a href="#">NM_001754.2</a></u> , <u><a href="#">NM_001754.3</a></u> , <u><a href="#">NM_001754.4</a></u> , <u><a href="#">BC050363</a></u> , <u><a href="#">BC110828</a></u> , <u><a href="#">BC136380</a></u> , <u><a href="#">BC136381</a></u> , <u><a href="#">BC144053</a></u> , <u><a href="#">BM149149</a></u> , <u><a href="#">BX475024</a></u> , <u><a href="#">NM_001754.5</a></u> |
| <b>UniProt ID:</b>       | <u><a href="#">Q01196</a></u>   |
| <b>Synonyms:</b>         | AML1; AML1-EVI-1; AMLCR1; CBF2alpha; CBFA2; EVI-1; PEBP2aB; PEBP2alpha  |
| <b>Components:</b>       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)  |
| <b>Quality Control:</b>  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.  |
| <b>Storage:</b>          | Store at -20°C.   |
| <b>Stability:</b>        | The primer mix is stable for one year from date of shipping.  |



[View online »](#)