

Product datasheet for **HP225989**

CD32 (FCGR2C) Human qPCR Primer Pair (NM_201563)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9103
Forward Sequence:	TCCAAGCCTGTGACCATCACTG
Reverse Sequence:	CCACTACAGCAGCAACAATGGC
ACCN:	<u>NM_201563</u> , <u>NM_201563.1</u> , <u>NM_201563.2</u> , <u>NM_201563.3</u> , <u>NM_201563.4</u> , <u>NM_201563.5</u> , <u>BC137397</u> , <u>BC137423</u> , <u>BP320976</u>
UniProt ID:	<u>P31995</u>
Synonyms:	CD32; CD32C; CDW32; FCG2; FCRIIC; IGFR2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)