

Product datasheet for **HP226010**

DPP10 Human qPCR Primer Pair (NM_020868)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57628
Forward Sequence:	GGACCGAAACCAAGTATGCTCTTC
Reverse Sequence:	CTGGAATCCAATTCCTCTGCCA
ACCN:	NM_020868 , NM_020868.1 , NM_020868.2 , NM_020868.3 , NM_020868.4 , BC030832 , BC030832.1 , BC038514 , BI597037 , BM974067
UniProt ID:	Q8N608
Synonyms:	DPL2; DPPY; DPRP-3; DPRP3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)