

Product datasheet for **HP226013**

Kallikrein 11 (KLK11) Human qPCR Primer Pair (NM_144947)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11012
Forward Sequence:	CGGCAACATCACAGACACCATG
Reverse Sequence:	CCAGGAGATAATGCCTTGAAGAG
ACCN:	<u>NM_144947</u> , <u>NM_144947.1</u> , <u>NM_144947.2</u> , <u>BC022068</u> , <u>BC022068.1</u> , <u>BC015551</u> , <u>BG697071</u> , <u>NM_144947.3</u>
UniProt ID:	<u>Q9UBX7</u>
Synonyms:	PRSS20; TLSP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)