

Product datasheet for **HP226035**

CDC25B Human qPCR Primer Pair (NM_021873)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	994
Forward Sequence:	AGAACCTCCTGGACAGTGACCA
Reverse Sequence:	GCTGAACTTGCCCGTCAATAGG
ACCN:	NM_021873 , NM_021873.1 , NM_021873.2 , NM_021873.3 , BC051711 , BC051711.1 , BC006395 , BC009953 , BM793120 , BX640836 , BX647988
UniProt ID:	P30305
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)